

Manual

Klenow fragment exo⁻

Fragment of DNA polymerase I from *E. coli*. Concentration 5 U/ μ l.

catalog #	size
1012-300	300 U

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

Description

Klenow fragment *exo*⁻ is the large fragment of DNA polymerase I from *E.coli*.

Enzyme possesses 5'-3' polymerase activity, but lacks 3'-5' exonuclease activity (proofreading) and 5'-3' exonuclease activity. 3'-5' exonuclease activity is eliminated by mutations in 3'-5' exonuclease active site.

Contents

	1012-300	storage
Klenow fragment <i>exo</i>⁻	300 U (5 U/ μ l)	-20 °C
storage buffer: 25 mM Tris-HCl, pH 8.0, 25 mM NaCl, 50% glycerol (v/v).		
Klenow reaction buffer	1 x 1.5 ml	-20 °C
10x PCR reaction buffer: 50 mM MgCl ₂ , 500 mM Tris-HCl, pH 8.0.		

Application

- random-primed DNA labeling
- labeling by fill-in 5'-overhangs of dsDNA
- strand displacement amplification (SDA)
- sequencing Sanger method
- pyrosequencing

Comments and recommendations

- incorporates modified nucleotides (e.g. Cyn3-, Cy5-, fluorescein-, rhodamine-, aminoallyl-, biotin-labeled nucleotides).
- inactivation by incubation for 10 min at 75 °C or by addition of EDTA.
- it is not recommended for DNA blunting reaction prior to DNA ligation since it frequently adds one or more extra nucleotides to 3'-terminus of blunt-end DNA substrates in non-templates directed fashion.
- the enzyme is active in restriction enzymes, PCR and RT buffers.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone: +48 883 323 761,+48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2023-2

