

Manual

RUN-HS 5 DNA polymerase

Hot Start *Taq* DNA polymerase with the reaction buffer. Anti-*Taq* monoclonal antibody-based Hot Start technology.
Concentration 5 U/ μ l.

catalog #	size
1001-200H-5	200 U
1001-1000H-5	1000 U

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

Advantages

- *Taq* DNA polymerase is the most popular DNA polymerase in PCR procedures.
- The enzyme is recommended for routine PCR reactions.
- Hot start technology avoids nonspecific amplification and enables room temperature reaction setup.

Description

RUN-HS 5 DNA polymerase is *Taq* polymerase purified from *E. coli* strain carrying a plasmid with a cloned gene encoding a DNA polymerase from *Thermus aquaticus*. Enzyme catalyzes incorporation of deoxynucleotides to 3' end of dsDNA at temperature 70-80 °C and presence of Mg^{2+} ions.

Polymerase is blocked with the anti-*Taq* monoclonal antibody. Full activation time requires 5 min of incubation at 95 °C. *Taq* DNA polymerase lacks 3'-5' exonuclease activity, but possesses weak 5'-3' exonuclease activity.

Contents

	1001-200H-5	1001-1000H-5	storage
RUN-HS 5 polymerase	200 U	1000 U	-20 °C

storage buffer:

100 mM KCl, 50 mM Tris-HCl pH 8.5, 0.5% Tween, 0.5% Triton X-100, 50% glycerol (v/v).

RUN-HS reaction buffer	1 x 1.5 ml	4 x 1.5 ml	-20 °C
------------------------	------------	------------	--------

10x PCR reaction buffer:

100 mM KCl, 100 mM $(NH_4)_2SO_4$, 200 mM Tris-HCl pH 8.5, 20 mM $MgSO_4$, 1% Igepal.

Unit definition

One unit of the enzyme catalyzes the incorporation of 15 nmol of dNTP into a polynucleotide fraction in 30 min at 75 °C.

Notes

- Before using, thoroughly thaw and gently mix by inverting the tubes.

PCR protocol

1. Add to the PCR tubes:

component	volume	final concentration
	50 μ l	
RUN-HS reaction buffer	5 μ l	1X
dNTP Mix (10 mM)	1-1.25 μ l	200-250 μ M
Primer 1 (10 μ M)*	1 μ l	0.2 μ M
Primer 2 (10 μ M)*	1 μ l	0.2 μ M
RUN-HS 5 polymerase	0.25 μ l	1.25 U
DNA template	variable	10 pg -1 μ g
ultrapure water	up to 50 μ l	

*For optimization, a primer titration should be performed from 0,2 μ M do 1 μ M final concentration.

2. Gently mix the samples and briefly centrifuge.

3. Place the tubes in the thermocycler and start the PCR programme.
An example amplification profile:

reaction step	temperature	time	number of cycles
enzyme activation	95 °C	5 min	1
denaturation	95 °C	15 s	
annealing*	50-68 °C	30 s	40
extension**	72 °C	30 s	

*Annealing temperature depends on primer sequence and the composition of the reaction mixture.

**Time of extension depends on the length of the amplicon.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone: +48 883 323 761,+48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2024-1

