

## *Manual*

# Plasmid Mini AX

Increased efficiency kit for low- and high-copy plasmid DNA purification.  
Procedure with DNA precipitation.

catalog #	size
010-50	50 isolations

For research use only.

### **Guarantee**

A&A Biotechnology provides a guarantee on this product.

The company does not guarantee the correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

# Table of Contents

<b>Specification</b>	<b>3</b>
<b>Contents</b>	<b>3</b>
<b>Additional equipment and reagents</b>	<b>3</b>
Necessary	3
Optional	3
<b>Important notes</b>	<b>4</b>
<b>Protocol</b>	<b>4</b>
<b>LySee color system</b>	<b>6</b>
Resuspension and lysis	6
Neutralization and precipitation	6
<b>Safety Information</b>	<b>7</b>

## Specification

<b>form</b>	midicolumn
<b>binding capacity</b>	20 µg of DNA
<b>sample size</b>	up to 10 ml of bacteria culture
<b>elution volume</b>	precipitation

## Contents

<b>component</b>	<b>size</b>	<b>storage</b>
Plasmid 20 columns	50 pcs	2–8 °C
20 ml tubes	50 pcs	15–25 °C
2 ml tubes	50 pcs	15–25 °C
L1 cell suspension solution	33 ml	2–8 °C
L2 lysis solution	33 ml	15–25 °C
L3T neutralizing solution	33 ml	15–25 °C
K1 equilibrating solution	55 ml	15–25 °C
K2P wash solution	220 ml	15–25 °C
K3 elution solution	70 ml	15–25 °C
PM precipitation mix	45 ml	15–25 °C
TE buffer	10 ml	15–25 °C

## Additional equipment and reagents

### Necessary

- 70% ethanol
- Centrifuge
- 2 ml sterile Eppendorf tubes

### Optional

- Sterile water (nuclease free) (cat.# 003-075, 003-25)

## Important notes

- Kit contains the LySee color system for easy optical control of alkaline lysis progress (page 6).
- SDS detergent is a component of L2 lysis solution and precipitates at low temperatures. Whenever the L2 lysis solution is not clearly transparent it must be warmed at 40 °C to form a thoroughly clear solution.

## Protocol

1. Centrifuge up to **10 ml** of overnight bacterial culture.

2. Discard the supernatant. Suspend the bacterial pellet in **600 µl** of **L1** cell suspension solution.

**Note.** During the pellet bacterial suspension, the solution will change color from a transparent deep pink to opaque light pink. The suspension is completed with complete disappearance of the pellet at the bottom tube.

Transfer the contents to 2 ml tubes (not included).

3. Add **600 µl** of **L2** lysis solution and gently mix. Keep for **3 min** at **room temp.**

**Note.** After the addition of L2 lysis solution, gently mix the tube so as not to cause fragmentation of the chromosomal DNA. Gently mix the tube by inverting a few times. The mixture should change appearance and color. After 3 min of incubation, the lysate must be completely clear and uniformly raspberry. If not, mix the lysate a few times and incubate again for 3 min at room temp.

4. Add **600 µl** of **L3T** neutralizing solution and gently mix until the disappearance of the raspberry color of the lysate.

**Note.** After the addition of L3T neutralizing solution followed by the rapid precipitation of the potassium salts (SDS), chromosomal DNA and certain proteins. After mixing, the tube contents should change the color to yellowish. No traces of raspberry color indicates complete neutralization and successful ending of the alkaline lysis.

5. Centrifuge for **5 min** at **10 000-15 000 RPM (~12 000 x g)**.

6. Place the Plasmid 20 column into a **20 ml tube** (included). Set the columns with tubes in the rack.

7. Apply **1 ml** of **K1** equilibrating solution onto the Plasmid 20 column. Wait for the solution to flow through the column.

8. Apply the clear lysate (supernatant) onto the Plasmid 20 column. Wait for the lysate to flow through the column.

9. Add **4 ml** of **K2P** wash solution. Wait for the solution to flow through the column.

10. Apply **200 µl** of **K3** elution solution directly to the Plasmid 20 column membrane. Wait for the eluate to flow through the column.

**Note.** The purpose of this step is to decrease the total volume of eluate, since the column void volume is about 200 µl.

11. Transfer the Plasmid 20 column to a **new 2 ml precipitation tube** ( included).

**Note:** The Plasmid 20 column drop director possesses proper fitting that allows easy attachment to the precipitation tube.

12. Add **1 ml of K3 elution solution**. Wait for the eluate to flow through the column.  
Remove the Plasmid 20 column.

13. PM precipitation mix contains a precipitation enhancer and it should be intensively mixed before use by vigorous hand shaking.

Add **800 µl of PM precipitation mix** to the eluted DNA.

14. Mix the sample by inverting the tube a few times and centrifuge for **10 min at 10 000 RPM (~10 000 x g)**.

15. Carefully discard supernatant. Be careful not to remove the DNA pellet at the bottom of the tube.

**Attention.** When pouring out the supernatant, it is very easy to lose the DNA pellet. For safety, it is recommended to pour the supernatant into the prepared tube so the pellet can be recovered.

16. Add **500 µl of 70% ethanol** (not included). Mix the sample and centrifuge for **5 min at 10 000 RPM (~10 000 x g)**.

**Note.** The light-blue DNA pellet should be visible at the bottom of the precipitation tube.

17. Carefully discard supernatant. Be careful not to remove the DNA pellet at the bottom of the tube.

**Attention.** When pouring out the supernatant, it is very easy to lose the DNA pellet. For safety, it is recommended to pour the supernatant into the prepared tube so the pellet can be recovered.

18. Air dry the plasmid DNA pellet for **5 min at room temp.** up-site down.

**Note.** If there are any leftovers (small droplets) of alcohol on the tube walls they should be removed with sterile cotton buds.

19. Dried DNA pellets can be dissolved in **50-150 µl of TE buffer** (included) or sterile water (not included).

**Note.** The blue color of DNA precipitate enables visual confirmation of the DNA dissolution process.

20. Store the plasmid DNA at **4-8 °C**.

## LySee color system

The LySee color system enables an easy and convenient visual control of alkaline lysis. The visual control system prevents common handling errors of incomplete cell resuspension, inefficient cell lysis and incomplete precipitation of unwanted cell components.

### Resuspension and lysis

The addition of the transparent purple L1 color cell suspension solution to the bacterial cell pellet makes the bacterial cell pellet easy to localize (fig 1). During the suspension of the bacterial cell pellet, the solution turns opaque light pink (fig 2). The suspension is completed with the complete disappearance of the pellet at the bottom of the tube. After the addition of L2 lysis solution and incubation, lysate turns transparent raspberry. Cell lysis is completed when the solution will turn homogeneously transparent raspberry (fig 3).




fig 1




fig 2



fig 3

### Neutralization and precipitation

The addition of the L3 neutralizing solution causes rapid precipitation of potassium salts (SDS), chromosomal DNA and some proteins (fig 4). After mixing, the solution turns yellowish (fig 5). No traces of raspberry color indicates complete neutralization and successful ending of alkaline lysis (fig 6).



fig 4



fig 5



fig 6

# Safety Information



**WARNING**

## L2 lysis solution

H315 Causes skin irritation.  
 H319 Causes serious eye irritation.  
 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
 H335 May cause respiratory irritation.  
 P261 Avoid breathing dust.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P342+P311 If experiencing respiratory symptoms call a Poison Center or doctor/physician.



**DANGER**

## L3T neutralizing solution

H315 Causes skin irritation.  
 H318 Causes serious eye damage.  
 H335 May cause respiratory irritation.  
 P261 Avoid breathing vapors.  
 P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



**WARNING**

## K1 equilibrating solution

H302 Harmful if swallowed.  
 H315 Causes skin irritation.  
 H319 Causes serious eye irritation.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



**DANGER**

## K3 elution solution

H225 Highly flammable liquid and vapor.  
 H319 Causes serious eye irritation.  
 H336 May cause drowsiness or dizziness.  
 P210 Keep away from heat, sparks, open flames, hot surfaces. No smoking.  
 P261 Avoid breathing vapors.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



**DANGER**

## Isopropanol

H225 Highly flammable liquid and vapor.  
 H319 Causes serious eye irritation.  
 H336 May cause drowsiness or dizziness.  
 P210 Keep away from heat, sparks, open flames, hot surfaces. No smoking.  
 P261 Avoid breathing vapors.  
 P305+P351+P338 If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



**A&A BIOTECHNOLOGY**  
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland  
phone +48 883 323 761, +48 600 776 268  
info@aabiot.com, www.aabiot.com

