

## *Manual*

# RUN-HS DNA polymerase

Taq DNA polymerase with reaction buffer.

Polymerase is blocked with anti-Taq monoclonal antibody (mAb). Form: solution. Concentration 1 U/ $\mu$ l.

catalog #	size
1001-200H	200 U
1001-1000H	1000 U

For research use only.

### **Guarantee**

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

## Description

**RUN-HS DNA polymerase** is *Taq* polymerase purified from *E.coli* stream carrying a plasmid with a cloned gene encoding a DNA polymerase from *Thermus aquaticus*.

Enzyme catalysis incorporation of deoxynucleotides to 3' end of dsDNA at temperature 70-80 °C and presence of Mg<sup>2+</sup> ions.

Polymerase is blocked with anti-Taq monoclonal antibody. Full activation time requires 3-5 min of incubation at 95 °C.

*Taq* DNA polymerase lacks 3'-5' exonuclease activity, but possesses weak 5'-3' exonuclease activity.

## Contents

	1001-200H	1001-1000H	storage
<b>RUN-HS polymerase</b>	200 U (1 U/μl)	1000 U (1 U/μl)	-20 °C
<b>RUN-HS reaction buffer</b>	1 x 1.5 ml	4 x 1.5 ml	-20 °C

10x PCR reaction buffer:  
100 mM KCl, 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 200 mM Tris-HCl, pH 8.5, 20 mM MgSO<sub>4</sub>, 1% Triton X-100.

## Notes

- Before using, thoroughly thaw and gently mix by inverting the tubes.

## Example PCR protocol

1. Thaw all components on ice, gently mix by inverting the tubes and briefly centrifuge. Place the tubes on ice again.
2. Place PCR tubes on ice and add:

component	PCR reaction volume
	50 $\mu$ l
RUN-HS reaction buffer	5 $\mu$ l
dNTP Mix (10 mM)	200-250 $\mu$ M (1-1.25 $\mu$ l)
Starter 1	0,1-0,5 $\mu$ M
Starter 2	0,1-0,5 $\mu$ M
RUN-HS polymerase	2-5 U
DNA template	10 pg - 1 $\mu$ g
Sterile water	up to 50 $\mu$ l

3. Gently mix the samples and briefly centrifuge.
4. Place the tubes in the thermocycler and start the PCR programme.

An example amplification profile for products up to 1000 bp:

step	temperature	time
Initial denaturation	95 °C	3-5 min
25-45 cycles	95 °C	15 s
	50-68 °C	30-60 s
	72 °C	1 min
Final incubation	72 °C	5-10 min

5. PCR products store in a refrigerator or freezer until later use.



**A&A BIOTECHNOLOGY**  
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland  
phone: +48 883 323 761,+48 600 776 268  
info@aabiotech.com, www.aabiotech.com

version 2023-2

