

Manual

E.coli Transformer Kit

Kit for preparation of competent *E.coli* cells and transformation. Chemical method.

| cat# | size |
|----------|------------------------|
| 4020-240 | 6 x 40 transformations |

For research use only.

Guarantee

A&A Biotechnology provides a guarantee on this product.

The company does not guarantee the correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components



Table of contents

| | |
|--|----------|
| Contents | 3 |
| Additional equipment and reagents | 3 |
| Necessary for competent cells preparation | 3 |
| Necessary for competent cells transformation | 3 |
| Media preparation | 4 |
| LB Agar (LA) (cat.# 2021-250, 2021-1000) | 4 |
| LB Miller (LB) (cat.# 2020-250, 2020-1000) | 4 |
| Competent cells preparation protocol | 5 |
| Competent cells transformation protocol | 6 |

Contents

| | 4020-240 | storage |
|--------------|----------|---------|
| S1E solution | 15 ml | +4 °C |
| S2E solution | 15 ml | +4 °C |

This kit was tested on derivatives of *Escherichia coli*: B, K-12 strains..

Additional equipment and reagents

Necessary for competent cells preparation

- *E.coli* strain
- sterile LB Miller medium (LB) (cat.# 2020-250, 2020-1000)
- sterile LB Agar medium (LA) (cat.# 2021-250, 2021-1000)
- sterile 1.5 ml Eppendorf tubes, sterile 50 ml Falcon tubes
- shaking incubator set to 37 °C
- centrifuge with rotor for 50 ml tubes

Necessary for competent cells transformation

- sterile LB Miller medium (LB) (cat.# 2020-250, 2020-1000)
- selection medium plates: 1 plate for 1 transformation
- sterile 50 ml Falcon tubes
- thermoblock set to 37 °C
- centrifuge with rotor for 50 ml tubes

Media preparation

Preparation of 1000 ml of medium.

LB Agar (LA) (cat.# 2021-250, 2021-1000)

1. Add **40 g** of medium to the appropriate vessel.

2. Add sterile water up to **1000 ml** and mix.

3. Autoclave for **10-20 min** at **121 °C**.

4. After cooling to **50-60 °C** mix again before use.

Note: At 25 °C pH should be 7.0.

LB Miller (LB) (cat.# 2020-250, 2020-1000)

1. Add **25 g** of medium to the appropriate vessel.

2. Add sterile water up to **1000 ml** and mix.

3. Autoclave for **10-20 min** at **121 °C**.

4. After cooling to **50-60 °C** mix again before use.

Note: At 25 °C pH should be 7.0.

Competent cells preparation protocol

- *Escherichia coli* should be streaked reductively onto LA medium
- Incubate plates overnight at 37 °C
- S1E and S2E solutions must be cooled on ice

1. Inoculate a single colony of *E.coli* obtained from reduction culture into **10 ml of LB medium** (if necessary with appropriate antibiotic).

Incubate **overnight at 37 °C**.
2. Add **1 ml of overnight culture** to **100 ml** of fresh **LB medium** (if necessary with appropriate antibiotic).
3. Incubate in a shaking incubator for **2-3 hours** at **220-250 RPM** at **37 °C** until $OD_{600}=0.3-0.5$ (logarithmic growth phase).
4. Keep the culture **on ice** for **15 min**.
5. Centrifuge for **5 min** at **3500 RPM** (~1500 x g) at **+4 °C**. Discard the supernatant.

If supernatant is not clear, increase centrifugation up to 3000 x g up to 10 min.
6. Resuspend the pellet in **2 ml** of **S1E solution**.
7. Carefully add **2 ml** of **S2E solution** and gently mix.
8. Keep **on ice** for **15 min**.
9. Transfer **100 µl** of **competent cell suspension** to 1.5 ml tubes.

100 µl of competent cells should be used for one transformation. Avoid multiple freeze-thaw cycles.
10. Competent cells are ready for transformation (page 5.) or can be stored at **-80 °C** for later use.

It is very important that competent cells are slowly frozen. Cells must not be frozen in liquid nitrogen.

Competent cells transformation protocol

- Prepare plates with a selective medium.
- Before transformation, the medium should be at room temp.
- It is recommended to prepare an additional plate with a selective medium for the negative control.

1. Use **100 µl** of ***E.coli* competent cells** for each transformation. Cells should be thawed on ice or freshly prepared and cooled on ice.
2. Add **plasmid DNA** or **ligation mixture** to the competent cells and gently mix.
The volume of DNA / ligation mixture should not exceed 20 µl.
3. Keep on ice for **15-45 min**.
During incubation avoid shaking or moving the tube. This reduces the effectiveness of the transformation. Longer incubation time improves the efficiency of the transformation.
4. Perform a heat shocking **60 s** at **42 °C**.
It is also possible to perform this step in an ultrasonic bath. Mixture should be placed in an activated ultrasonic bath for 5-10 s at room temp.
5. Keep on ice for **2 min**.
If transformed plasmid is ampicillin resistant follow point 7.
6. If transformed plasmid is not ampicillin resistant (or with resistance for another antibiotic):
 - add **1 ml** of **preheated to 37 °C LB media without antibiotics**.
 - incubate for **45 min** with shaking at **220 RPM** at **37 °C** to express genes responsible for antibiotic resistance.
7. Cultivate **20-200 µl** of **transformation mixture** on a plate with a selective medium.
8. Incubate **overnight** at **37 °C**.

Average efficiency of transformation of *E.coli* TOP10F' cells with pUC19 plasmid is 10⁷-10⁸ CFU.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone +48 883 323 761, +48 600 776 268
info@aabiotech.com, www.aabiotech.com

