

Manual

RUN DNA polymerase

Taq DNA polymerase with reaction buffer. Concentration 1 U/ μ l.

catalog #	size
1001-200	200 U
1001-1000	1000 U

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components



Advantages

- *Taq* DNA polymerase is the most popular DNA polymerase in PCR procedures.
- Recommended for routine PCR reactions.

Description

RUN DNA polymerase is *Taq* polymerase purified from *E.coli* strain carrying a plasmid with a cloned gene encoding a DNA polymerase from *Thermus aquaticus*.

Enzyme catalysis incorporation of deoxynucleotides to 3' end of dsDNA at temperature 70-80 °C and presence of Mg²⁺ ions.

Taq DNA polymerase lacks 3'-5' exonuclease activity, but possesses weak 5'-3' exonuclease activity.

Contents

	1001-200	1001-1000	storage
RUN polymerase	200 U (1 U/μl)	1000 U (1 U/μl)	-20 °C
storage buffer: 10 mM KCl, 20 mM Tris-HCl pH 8.7, 0,1 mM EDTA, stabilizers, 50% glycerol (v/v).			
RUN reaction buffer	1 x 1.5 ml	4 x 1.5 ml	-20 °C
10x PCR reaction buffer: 100 mM KCl, 100 mM (NH ₄) ₂ SO ₄ , 200 mM Tris-HCl, pH 8.5, 20 mM MgSO ₄ , 1% Triton X-100.			

Notes

- Before using, thoroughly thaw and gently mix by inverting the tubes.

Example PCR protocol

1. Thaw all components on ice, gently mix by inverting the tubes and briefly centrifuge. Place the tubes on ice again.
2. Place PCR tubes on ice and add:

component	PCR reaction volume
	50 μ l
RUN reaction buffer	5 μ l
dNTP Mix (10 mM)	200-250 μ M (1-1.25 μ l)
Starter 1	0,1-0,5 μ M
Starter 2	0,1-0,5 μ M
RUN polymerase	1-2 U
DNA template	10 pg - 1 μ g
Sterile water	up to 50 μ l

3. Gently mix the samples and briefly centrifuge.
4. Place the tubes in the thermocycler and start the PCR programme.

An example amplification profile for products up to 1000 bp:

step	temperature	time
Initial denaturation	94 °C	1-5 min
25-45 cycles	94 °C	30-60 s
	50-68 °C	30-60 s
	72 °C	1 min
Final incubation	72 °C	5-10 min

5. PCR products store in a refrigerator or freezer until later use.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone: +48 883 323 761,+48 600 776 268
info@aabiotech.com, www.aabiotech.com

version 2023-2

