

Manual

dNTP Mix 10 / dNTP Mix 25

Mix of deoxynucleotides (dATP, dCTP, dGTP, dTTP) at concentration of 10 mM / 25 mM of each dNTP.

catalog #	size
2001-200	200 μ l
2001-1000	1 ml
2001-1	1 ml
2001-4	4 x 1 ml

For research use only.

Guarantee

A&A Biotechnology provides guarantee on this product.

The company does not guarantee correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components

Description

dNTP mix contains a tube with a water solution of dATP, dCTP, dGTP, dTTP) at a concentration of:

10 mM # 2001-200, # 2001-1000

25 mM # 2001-1, # 2001-4

Highest purity.

Store at -20 °C

Thaw and carefully mix before use.



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, ul. Strzelca 40, 80-299 Gdańsk, Poland
phone +48 600 776 268, +48 883 323 761
info@aabiotech.com, www.aabiotech.com

