



A&A BIOTECHNOLOGY
innovating life science

Manual

Ultrapure water

Molecular biology grade. Nuclease free. Not treated with DEPC.

catalog #	size	catalog #	size
005-15	1.5 ml	005-5	5 ml
005-515	5 x 1.5 ml	005-55	5 x 5 ml
005-1015	10 x 1.5 ml	005-100	100 ml
005-2515	25 x 1.5 ml	005-250	250 ml
005-10015	100 x 1.5 ml	005-500	500 ml

For research use only.

Guarantee

A&A Biotechnology provides a guarantee on this product.

The company does not guarantee the correct performance of this kit in the event of:

- not adhering to the supplied protocol
- use of not recommended equipment or materials
- use of other reagents than recommended or which are not a component of the product
- use of expired or improperly stored product or its components



Table of Contents

Description	3
Specification	3

Description

Ultrapure water is ideal for use in all molecular biology applications.

It is rigorously tested for endonucleases, exonucleases, DNAses and RNAses contamination.

Specification

- Not DEPC-treated
- Reverse osmosis filtered
- Deionized
- UV treated
- Membrane filtered



A&A BIOTECHNOLOGY
innovating life science

A&A Biotechnology, al. Zwycięstwa 96/98, 81-451 Gdynia, Poland
phone +48 58 6982194, fax +48 58 6228578
info@aabiotech.com, www.aabiotech.com

